本项目地址: https://github.com/soniccomm/sonic_application

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https://github.com/soniccomm/sonic_application

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Overview

This repository provides a comprehensive MATLAB-based toolkit for ultrasound imaging processing, specifically designed to handle raw data acquired via **plane-wave ultrasound** transmission. It includes post-processing scripts for four advanced blood flow imaging modalities, enabling the generation of rich hemodynamic images from raw data.

Core Functional Modules

- **UltrafastColorDoppler**: Generates color flow maps depicting blood velocity and direction.
- **UltrafastPowerDoppler**: Produces power Doppler images highlighting the presence and relative density of blood flow, with enhanced sensitivity to low-velocity flow.
- **UltrafastSpectralDoppler**: Extracts velocity spectra at specific points for quantitative hemodynamic analysis.
- **UltrasoundLocalizationMicroscopy (ULM)**: Achieves super-resolution microvascular imaging beyond the diffraction limit by tracking ultrasound contrast agents (microbubbles).

Preparations: Data Acquisition

Critical Note: All scripts are **post-processing tools**. Users **must** have already **acquired and saved raw plane-wave data** before running these scripts.

Download & Setup

1. Download the Repository:

- Clone or download the entire project repository from GitHub
Link: https://github.com/soniccomm/sonic_application

2. Folder Structure Setup:

- Place the Application folder into the script directory.
- Move all files from the lib folder into the lib subfolder within the script directory.

3. **MATLAB Configuration:**

- Navigate MATLAB's *Current Folder* to the specific script directory you intend to use (e.g., .../Application/UltrafastPowerDoppler/).
- Open the target script and modify the user-defined parameters in the configuration section (e.g., data paths, acquisition settings) to match your dataset.
- Run the script to execute the processing pipeline.